## Storage Requirements

Store between 5° C and 25° C.

## Product Code

OP8G

## Intended Use

ODONTOPASTE® is used to preserve the disinfection of a root canal. This disinfection may be achieved with sodium hypochlorite and EDTA rinses along with standard endodontic filing.

ODONTOPASTE® is applied following the disinfection procedures. ODONTOPASTE® acts as an obturator of the root canal by occupying any space which may allow for bacteria to harbour or repopulate. In doing so, ODONTOPASTE® allows for healing and hence a reduction in inflammatory exudates so that a dry surface can be presented to the adhesives and sealers that are used to seal the canal. ODONTOPASTE® is intended for continuous use for less than 30 days.

## Instructions for Use

ODONTOPASTE® is a water soluble, white, non-staining paste. ODONTOPASTE® contains clindamycin hydrochloride and triamcinolone acetonide. Clindamycin hydrochloride is a broad-spectrum bacteriostatic antibiotic while triamcinolone acetonide is a steroidal anti-inflammatory agent. The clindamycin hydrochloride acts as a preservative for the bulk of the obturating paste, preventing bacterial repopulation within the paste. The triamcinolone acetonide provides temporary pain relief.

Only when used in conjunction with proper endodontic procedures will ODONTOPASTE® assist with the relief of pain associated with acute pulpitis and periodontitis. ODONTOPASTE® is used as an interim pulpotomy dressing material for permanent teeth in cases of acute pulpitis where conventional pulpectomy and root canal treatment will follow.

ODONTOPASTE® cannot kill bacteria as the antibiotic is bacteriostatic and NOT bactericidal. Any longstanding infections will be unaffected by ODONTOPASTE® unless the systemic immune system has the capacity to infiltrate the biological cause. This is difficult in avascular root canals. For such infections, a calcium hydroxide-based root canal paste is recommended for its active antibacterial action.

ODONTOPASTE® is designed to be used as a therapeutic dressing material in prepared canals after pulpectomy and root canal preparation to

- temporarily reduce the pain of apical periodontitis
- reduce bacterial repopulation of root canals

ODONTOPASTE® can also be used in prepared root canals following trauma to reduce the incidence of inflammatory root resorption and in the reduction of existing inflammatory root resorption.

ODONTOPASTE® has negligible to no intrinsic staining of teeth, is not photosensitive and, as such, is not limited to subgingival application.

The device can be used on all age groups and body types.

ODONTOPASTE® should only be applied by dentists.

## Method of Use

- Prepare the tooth cavity according to the standard dental procedure for root canal medication application. Ensure the ODONTOPASTE® has not exceeded its expiration date.
- Remove from the packaging, ensure the tamper evident caps are still intact before its first use, unscrew cap.
- Dispense the required amount (sufficient quantity to fill the tooth cavity) onto a dental mixing pad.
- Screw up cap, record when the tube was opened.
- Apply the ODONTOPASTE® directly onto inflamed pulp with a cotton pellet and/or a lentulo spiral (dental tool)
- Ensure the root canal cavity is filled with the ODONTOPASTE®.
- Seal the tooth as per standard dental procedure for applying a dental restoration.
- Ensure the entire tube of ODONTOPASTE® is utilised prior to the expiration date.
- More than one application may be required; the need for re-application should be decided by a dental professional.

## Contra-indications for Use

- The ODONTOPASTE® is supplied non-sterile; the device is not designed to be sterile.
- Each tube is designed to be used for numerous applications.
- The ODONTOPASTE® should not be used if it has exceeded its expiration date.
- Mixing surfaces and dental instruments should be cleaned in accordance with cross-infection protocols. Refer to relevant infection control protocols for dentistry.
- ODONTOPASTE® should not be mixed with additional calcium hydroxide as the steroid and antibiotic activity is impaired.

## Warnings

- Each user should read the instructions for use and seek assistance with initial use if required.
- Safety glasses, mask and gloves should be worn when handling the ODONTOPASTE®.
- ODONTOPASTE® is not recommended for use in patients with allergy to clindamycin or lincomycin.
- ODONTOPASTE® should not be used on patients taking concurrent doses of erythromycin as antagonism has been demonstrated between the two drugs in vitro.
- ODONTOPASTE® should be used with caution on patients with a history of gastrointestinal disease, particularly those with colitis.
- ODONTOPASTE® should be used with caution on nursing mothers as clindamycin has been shown to be present in breast milk.
- ODONTOPASTE® should be used with caution on pregnant mothers as triamcinolone acetonide has demonstrated teratogenic effects in pregnant animals.

## Disposal

- Please dispose of unused portions in accordance with all applicable local, state and federal regulations.
- Empty tube is recyclable aluminium.
- Nozzle should be discarded after each tube.

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